

OPTIMA MODULAR FEMORAL COMPONENT *for Hemi-Arthroplasty*

ENDO HEADS

– Nine diameters in two neck lengths

PRECISE taper/trunnion machining

PLATFORM PERPENDICULAR
TO RESULTANT FORCES

– to maximize proximal loading

PLATFORM AND "COBRA" FLANGES

– to increase cement pressurization

LATERAL FIN

– to enhance rotational stability

COBALT-CHROMIUM STEM

DISTAL TAPER

– to reduce incidence of distal thigh pain
and endosteal wall pressure

SIMPLE INSTRUMENTATION

ECONOMICAL

NO POLYETHYLENE WEAR DEBRIS

BIOLOGICALLY ORIENTED PROSTHESES
BIOPRO



Optima Modular Femoral Component



Optima Modular Hip Stem

Size 2	10160
Size 4	10161
Size 6	10162
Size 8	10163

Endo Heads

Short Neck

Ceramic	Cobalt Chrome	Size
13006	10188	38
13007	10189	41
13008	10190	43
13009	10191	45
13010	10192	47
13011	10193	49
13012	10194	51
14823	10195	53
14825	10196	55

Medium Neck

Ceramic	Cobalt Chrome	Size
13013	10179	38
13014	10180	41
13015	10181	43
13016	10182	45
13017	10183	47
13018	10184	49
13019	10185	51
14824	10186	53
14826	10187	55

*For more information contact your representative
or call Biopro at 1-800-252-7707*



www.biopro-imt.com